# CONAN

## Instructions for use

### Item

Glove nitrile powder-free  
Conan 10-019  
10x100 pieces  
Blue | Green | Black  
6/XS | 7/S | 8/M | 9/L | 10/XL | 11/2XL

### Standards

Category III

EN ISO 374-1:2016+A1:2018  
Type B -KPT

EN ISO 374-5:2016  
VIRUS

EN 455

(EU) No 10/2011

(EC) No 1935/2004

Notified body responsible for certification and ongoing conformity:

SATRA Technology Europe Limited  
Bracetown Business Park, Clonee, D15YN2P, Republic of Ireland

European authorized representative:

Lotus NL B.V.  
Koningin Julianaplein 10, 1<sup>st</sup> floor, 2595 AA, The Hague, The Netherlands, peter@lotusnl.com

### SPC B.V.

PO Box 15700  
1001 NE Amsterdam  
The Netherlands

Available in the following languages below, click on the box to select your preferred language:

EN

DE

NL

FR

PL

| Version | Date     | Revision date |
|---------|----------|---------------|
| 1.0     | 1-9-2025 | 1-9-2025      |

## Instruction for use

EN

### Instruction for use for disposable examination and protective gloves made of nitrile (powder-free, non-sterile)

**Usage:** only wear the products in a suitable size. The optimal level of protection will not be provided if the glove is too loose or too tight. The gloves shall not be worn when there is a risk of entanglement with moving parts of machines. Before usage, inspect the gloves for any defect or imperfections. Do not use damaged gloves. People with allergies such as a skin rash or abnormal feeling, should immediately stop wearing the glove. It is the employer's responsibility together with the user to analyze if each glove protects against the risks that can appear in any given work situation. Please read the instructions for use, prior to using the products.

**CE:** the glove is in conformity with the provisions of Regulation (EU) 2016/425 (Category III) and the standard EN ISO 21420:202.

**Declaration of Conformity:** can be found on [www.conanprotection.com](http://www.conanprotection.com).

**The gloves are designed to protect against the following risk:**

#### EN ISO 374-1:2016+A1:2018 - Protective gloves against dangerous chemicals and micro-organisms

**Penetration:** the glove shall not leak water or air when tested according to penetration, EN ISO 374-2:2019.

**Degradation:** indicates the change in puncture resistance after exposure to the challenge chemical. Degradation shall be determined according to EN ISO 374-4:2019 for each chemical.

**Permeation:** the glove must withstand a breakthrough time of at least:

Type A - 30 minutes (level 2) against minimum 6 test chemicals

Type B - 30 minutes (level 2) against minimum 3 test chemicals

Type C - 10 minutes (level 1) against minimum 1 test chemical

The test chemicals are listed in the table below and all 18 chemicals shall be tested according to EN 16523-1:2015+A1:2018.

**Micro-organisms:** the glove is tested to protect against bacteria, fungi and, if applicable, viruses, EN ISO 374-5:2016.

**Warning:** chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based from samples taken 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

**For single use only**

#### EN ISO 374-1: 2016+A1:2018

| Test chemical           | EN ISO 374-1:2016+A1:20018/Type B<br>Permeation level | EN ISO 374-4:2019<br>Degradation (mean value) |
|-------------------------|---|---|
| K Sodium Hydroxide 40%  | 6   | -52.7%  |
| P Hydrogen Peroxide 30% | 2   | 19.4%   |
| T Formaldehyde 37%      | 4   | 1.4%  |

| Performance level              | 1   | 2   | 3   | 4    | 5    | 6    |
|--------------------------------|-----|-----|-----|------|------|------|
| Min. breakthrough times (mins) | >10 | >30 | >60 | >120 | >240 | >480 |

#### EN ISO 374-5:2016

Resistance to bacteria and fungi: pass

Resistance to virus: pass

The penetration resistance had been assessed under laboratory conditions and relates only to the tested specimen.

## Instruction for use

### EN

#### Food Contact Materials Regulation

Suitable for contact with foodstuffs. Latex free, powder-free and DEHP/DOP-free. The glove is in conformity with the European Regulations (EU) No 10/2011 and (EC) No 1935/2004.

#### Medical Device Regulation

Compliant with MD regulations (EU) 2017/745 (Class I), EN 455-1:2020, EN 455-2:2015, EN 455-3:2015, EN 455-4 :2009.

#### Storage conditions

Store the gloves in their original packaging in a dark, cool and dry place with relative humidity  $\leq 80\%$  and temperature below  $40^{\circ}\text{C}$ , free of corrosive gases and well vented. The mechanical properties of the glove will not be affected when stored properly.

#### Proper donning & doffing

Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers and others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. **How to don gloves:** remove one glove from the package and inspect it to be sure no pinholes or tears are present. If gloves are ambidextrous, they can be worn on either hand. If not, align the glove's fingers and thumb with the proper hand before donning. Insert five fingers into the cuff and pull the cuff over the wrist. Check for a secure fit around the fingers and palm. The cuff should fit snugly around the wrist. **How to doff gloves:** when removing the gloves, avoid allowing the outer surface of the gloves to come in contact with your skin, because the surface may have been contaminated with blood and other body fluids. Avoid snapping, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby. Grasp the outside edge of the glove near the wrist. Peel the glove away from the hand, turning it inside out. Hold it in the opposite gloved hand. Slide an ungloved finger under the wrist of the remaining glove, being careful not to touch the outside of the glove. Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Discard.

**Cleaning:** the gloves are single-use only and not designed to be laundered nor to be reused. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization has not been validated for this product by the manufacturer. Product integrity may be compromised by any cleaning or re-sterilization process used.

**Disposal:** used products which have been in contact with chemicals or contaminated with infectious or other hazardous materials such as residual pesticides should be disposed after each working shift and not reused. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves. Dispose of according to Local Authority Regulations. Landfill or incinerate under controlled conditions. If the glove has been used in a clinical setting after use, the product should be incinerated or disposed as per the clinical waste disposal protocol.

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