

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Peracetic acid 50

Parameter*	Value*
REF	91340
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Peracetic acid
Gradation	0 • 5 • 10 • 20 • 30 • 50 mg/L Peracetic acid
Method	Redox reaction
Color change	white → blue
Instrumental measuring range	5–50 mg/L Peracetic acid

*These data are typical values and must not be regarded as specification.

